


# SWAB IDENTIFICATION CHART

## MICROBIOLOGY



- 1** Amies swab for bacterial culture

## CYTOLOGY



- 2** Surepath vial with broom and brush for cervical cytology  
(High risk HPV performed from vial as appropriate)



- 3** Copan 552C HPV Self Collect Swab

## STI NAAT

(Chlamydia, Gonorrhoea, Trichomonas and M. Genitalium ONLY)



- 4** Aptima Multi-Test and Unisex Swab

- 5** Aptima Urine Collection Kit

## VIRAL AND BORDETELLA PERTUSSIS PCR

(Including Herpes, Varicella, Measles and SARS-CoV-2)



- 6** Multiple Brands  
with Universal Transport Media (UTM) **and/or** Viral Transport Media (VTM)

## LABELLING

Ensure all specimens are labelled with the patient's first name, last name, NHI, date of birth and date/time of collection